

PKS™ PlasmaSORD

BIPOLAR MORCELLATOR

PK TECHNOLOGY®

The First Completely
Bladeless SORD
(SOLID ORGAN REMOVAL DEVICE)



- **Safe bipolar energy** to morcellate tissue instead of sharp, spinning blades.
- **Efficient bipolar energy** creates fewer tissue fragments to clean up.
- **Lightweight, ergonomic** handle reduces hand fatigue in prolonged cases.
- **Easy to use** "plug and play" PK® Technology with the G400 Generator requires little setup for O.R. staff.
- **Durable and economical** because only one bipolar device is required per case.

PKS™ PlasmaSORD™ Bipolar Morcellator

OLYMPUS

Your Vision, Our Future

ORDERING INFORMATION

PKS™ PLASMASORD

Order #	Description
962000PK	PKS™ PlasmaSORD™

PK * TECHNOLOGY G400 & ACCESSORIES

777000	G400 Generator w/9-Pin Inputs
744010	G400 Footswitch
560085-002	Generator Power Cord
700000	PKS™ and G400 Generator Cart

PK * TECHNOLOGY INSTRUMENTS

956000PC	PKS™ Plasma J-Hook,™ 15 cm
956005PC	PKS™ Plasma J-Hook,™ 24 cm
956010PC	PKS™ Plasma J-Hook,™ 33 cm
956015PC	PKS™ Plasma J-Hook,™ 45 cm
952000PK	PKS™ PlasmaSpatula,® 24 cm
952005PK	PKS™ PlasmaSpatula,® 33 cm
952010PK	PKS™ PlasmaSpatula,® 45 cm
921015PK	PKS™ Cutting Forceps, 15 cm 10mm

PK * TECHNOLOGY INSTRUMENTS, CONTINUED

920000PK	PKS™ Cutting Forceps, 24 cm 5mm
920005PK	PKS™ Cutting Forceps, 33 cm 5mm
921020PK	PKS™ Cutting Forceps, 33 cm 10mm
920010PK	PKS™ Cutting Forceps, 45 cm 5mm
950000PK	PKS™ Needle, 33 cm
942000PK	PKS™ LYONS™ Dissecting Forceps, 24 cm
942005PK	PKS™ LYONS™ Dissecting Forceps, 33 cm
942010PK	PKS™ LYONS™ Dissecting Forceps, 45 cm
954000PK	PKS™ L-Hook,™ 33 cm
954005PK	PKS™ L-Hook,™ 45 cm
910000PK	PKS™ MOLLY® Forceps, 33 cm 3mm
910005PK	PKS™ MOLLY® Forceps, 33 cm 5mm
910010PK	PKS™ MOLLY® Forceps 45 cm 5mm
960000PK	PKS™ LP Scissors, 33 cm
960005PK	PKS™ LP Scissors, 45 cm
915000PK	PKS™ SEAL® Open Forceps, Curved Jaw, w/Cord, 25 cm
915005PK	PKS™ SEAL® Open Forceps, Straight Jaw, w/Cord, 25 cm
915010PK	PKS™ SEAL® Open Forceps, Angled Jaw, w/Cord, 25 cm

SPECIFICATIONS

PKS™ PLASMASORD™

Shaft Length:	6.69 in. (170 mm)
Shaft Diameter:	0.59 in. (15 mm)
Grasper Diameter:	0.39 in. (10 mm)
Depth Gauge:	
Length	0.79 in. (20 mm)
Diameter	≤ 1.57 in. (40 mm)

PKS™ PLASMASORD™ CORD

Length:	12 ft. (3.65 m)
---------	-----------------

PK * G400 GENERATOR

Voltage:	Nominal 100–120 / 230 Volts RMS 50/60 Hz, 600 VA Source Operating Range
Waveform:	315–480 kHz
Dimensions (W x D x H):	16 x 16 x 5.3 in. (410 x 410 x 135 mm)
Weight:	17.6 lbs (8.0 kg)
Electro-medical:	Class I (IEC 60601-1), Class IIb (MDD)
Maximum output:	
PK & PC Modes	200 W into 150 ohms
DES Modes	128 W into 70 ohms
VP Modes	90 W into 80 ohms

PK * G400 FOOTSWITCH

Dimensions (W x H x D):	13.5 x 2.25 x 9 in.
Weight:	4.4 lbs (2.0 kg)
Classification:	IPX8 watertight type (except plug)

PK TECHNOLOGY including the PKS™ PlasmaSORD™ is
The Most Comprehensive Bipolar Energy Platform Available for LSH & TLH



PKS™ Cutting Forceps

The unique five-in-one design grasps, dissects, coagulates, transects and retracts tissue all in one easy-to-use device across a wide variety of procedures.

PKS™ Lyons™ Dissecting Forceps

Two-tiered, tapered jaw design with serrated surfaces provides secure grasping and retraction of tissue as well as pinpoint or wide based coagulation.

PKS™ Plasma Spatula®

Unique angled tip and broad surface provide the ability to precisely cut, coagulate and dissect in a variety of procedures.

PKS™ Plasma J-Hook

Simultaneously cuts and coagulates to reduce procedure time. Tip design is excellent for skeletonization & mobilization.

PKS™ Needle

Adjustable needle electrode allows the surgeon to adjust the depth of the cut for precise control and minimal damage to surrounding tissue.

For more information call 508.804.2600
For order inquiries call 888.524.7266
Visit us at www.gyrusacmi.com

SURGICAL AND ENDOSCOPY:

Gyrus ACMI
136 Turnpike Road
Southborough, MA 01772
TEL 508-804-2600
FAX 508-804-2624

Sales Locations:

USA

Gyrus ACMI
6655 Wedgwood Rd, Suite 160
Maple Grove, MN 55311
TEL 763-416-3000
FAX 763-463-1200

Canada

Olympus Canada, Inc.
151 Telson Road
Markham, Ontario L3R 1E7
Canada
TEL 1-800-387-0437
FAX 1-905-479-2595

UK

KeyMed Ltd
KeyMed House
Stock Road
Southend-on-Sea
Essex SS2 5QH
United Kingdom
TEL +44 (0)1702 616333
FAX +44 (0)1702 465677

Australia

Olympus Australia Pty Ltd.
31 Gilby Road
Mount Waverley
Victoria 3149
Australia
TEL 61 1300-132-992
FAX 61 (0) 3-9543-1350

Europe

Olympus Medical Systems Europa GmbH
Wendenstr. 14–18
20097 Hamburg
Germany
TEL + 49 40 23773-0
FAX + 49 40 23773-4656


GYRUS ACMI
Power through integration
– Olympus