



Better Grasping.
Superior Dissection.
Precise Coagulation.



LYONS™ DISSECTING FORCEPS
IS™ DISSECTING FORCEPS

GYRUS
ACMI

GRASP, DISSECT AND COAGULATE WITH ONE VERSATILE INSTRUMENT

The PKS LYONS™ Dissecting Forceps enables you to grasp, dissect and coagulate tissue quickly and effectively with a single instrument, reducing instrument exchanges and procedure time. With its unparalleled 3-in-1 versatility, the PKS LYONS™ Dissecting Forceps is the one instrument no laparoscopic surgeon should be without.

▲ BETTER GRASPING

Two-tier jaw design with serrated surfaces provides secure grasping and retraction of tissue.

▲ SUPERIOR DISSECTION

Finely rounded tip provides atraumatic dissection and assists with vessel skeletonization.

▲ PRECISE COAGULATION

Tapered jaw design combined with PK energy provides pinpoint or wide based coagulation while virtually eliminating sticking and thermal spread.

▲ PK TECHNOLOGY ADVANTAGE

One energy source powers multiple instruments and facilitates an entire procedure. Plug and play feature recognizes each instrument and optimizes output characteristics.

GRASP
DISSECT
COAGULATE



Pinpoint Coagulation



Wide Base Coagulation



ORDERING INFORMATION

DIAMETER	DESCRIPTION	LENGTH	ORDER #
<i>For use with Gyrus PK SuperPulse or PK Generator</i>			
5mm	PKS LYONS™ Dissecting Forceps	24cm	3720PK
5mm	PKS LYONS™ Dissecting Forceps	33cm	3700PK
5mm	PKS LYONS™ Dissecting Forceps	45cm	3740PK



Gyrus Medical, Inc.
6655 Wedgwood Road, Suite 160
Maple Grove, MN 55311-3602

For ordering information call: 800.852.9361
Outside US call: 763.416.3000
Fax: 763.463.1200

www.gyrusmedical.com
Email: info@gyrusmed.com

©2005 Gyru Medical, Inc. All rights reserved
6664 1005

Gyrus International, LTD.
410 Wharfedale Road, Winnersh Triangle
Wokingham, Berkshire RG41 5RA UK

Tel: 44 (0) 1189 219700
Fax: 44 (0) 1189 219800

www.gyrusinternational.com
Email: sales@gyrusintl.com

LYONS is a trademark of Gyru Medical, Inc.

